



Heat inactivated *Streptococcus oralis*

Cat. No.: LBGF-0224-GF48

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	<i>Streptococcus oralis</i> has been inactivated by heating to 65°C for 30 minutes.
Target	<i>Streptococcus</i> (Inactivated)
Product Type	Inactivated organism
Particle Integrity	Whole, intact particles
Risk Group	1
CFU	≥10 ⁶ CFU/mL; ≥10 ⁷ CFU/mL; ≥10 ⁸ CFU/mL
Volume	0.5 mL
Pack Size	1 vial
Product Format	Frozen
Storage	-20°C or colder
Shipping	Dry ice
Applications	<ol style="list-style-type: none">1. Extraction of nucleic acid reference materials;2. Specificity analysis of nucleic acid detection kit, cross reaction verification;3. Test method evaluation and laboratory quality control.
Alternative Names	Heat inactivated; Heat killed; <i>Streptococcus oralis</i>