



## Recombinant pgi protein from *Streptococcus pyogenes* M3

Cat. No.: LBGF-0626-G384

This product is for research use only and is not intended for diagnostic use.

### BASIC INFORMATION

#### Product Overview

Our Recombinant pgi protein is sourced from *Streptococcus pyogenes* M3 and expressed via specialized heterologous hosts (*E. coli*, yeast, or mammalian cells) to achieve native-like conformation. We prioritize exceptional purity and low endotoxin levels to guarantee reliable performance in your experimental workflows. Detailed quality metrics, including purity and specific biological activity (where available), can be found on the lot-specific Certificate of Analysis (CoA).

#### Target

pgi protein

#### Source Organism

*Streptococcus pyogenes* serotype M3

#### Gene/Protein Name

pgi

#### Protein Sequence

MSHITFDYSK VLESFAGQHE IDFLQGQVTE ADKLLREGTG PGSDFLGWLD LPENYDEE  
EF

ARILTAAEKI KSDSEVLVVI GIGGSYLGAK AAIDFLNHHF ANLQTAKERK APQILYAGNS

ISSTYLADLV EYVQDKEFSV NVISKSGTTT EPAIAFRVFK ELLVKKYGQE EASKRIYATT

DKVKGAVKVE ADANNWETFV VPDNVGGRFS VLTAVGLLPI AASGADITAL MEGANAAR  
KD

LSSDKISENI AYQYAAVRNL LYRKGYITEI LANYEPSLQY FGEWWKQLAG ESEKDKQKGI

YPTSANFSTD LHSLGQFIQE GYRNLFFETVI RVDNPRKNVI IPELAEDLDG LGYLQGKDV

D

FVNKKATDGV LLAHTDGGVP NMFVTLPAQD EFTLGTYIYF FELAIIVSGY MNAVNPFDDQ  
P

GVEAYKRNMF ALLGKPGFEA LSAELNARL

---

**Subcellular Localization** Cytoplasm

---

**Protein Region** Full Length

---

**UniProt Accession Number** P0DB24

---

**Biological Activity** Not Test

---

**Protein Family** GPI family

---

**Applications** ELISA, Western Blot (WB), Immunoprecipitation (IP), Positive Control.  
Note: Optimal dilutions/concentrations should be determined by the end user.

---

## CHEMICAL & PHYSICAL PROPERTIES

**Tag&Position** Determined during the production process by default. If a specific tag is required, please notify us, and its development will be prioritized.

---

**Purity** ≥85% (SDS-PAGE)

---

**Formulation**

1. Liquid Form: Supplied in a default Tris/PBS-based storage buffer supplemented with 5% -50% glycerol.
2. Lyophilized Form: Lyophilized from a Tris/PBS-based pre-lyophilization buffer (pH 8.0) containing 6% Trehalose as a lyoprotectant.

---

**Format** Liquid or Lyophilized powder

---

## HANDLING & STORAGE

**Stability&Storage**

1. Lyophilized Powder: Stored at -20°C to -80°C, the product is stable for 12 months from the date of receipt. Upon reconstitution, it is stable for up to 1 week at 2-8°C or 3 months at -20°C to -80°C.
2. Liquid Form: Stored at -20°C to -80°C, the product is stable for 6 to 12 months (avoid repeated freeze-thaw cycles). Aliquot upon arrival to prevent degradation.

---

|                       |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                       |
|-----------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Reconstitution</b> | <p>For Lyophilized Powder. Always centrifuge tubes before opening to bring the contents to the bottom. Do not vortex.</p> <ol style="list-style-type: none"> <li>1. Reconstitute in sterile ddH<sub>2</sub>O to a stock concentration of 0.1-1.0 mg/mL.</li> <li>2. Dissolve the protein by gentle pipetting up and down or by inverting the tube.</li> <li>3. Allow the solution to stand at room temperature for at least 20 minutes to ensure complete reconstitution.</li> <li>4. Aliquot the reconstituted solution into working volumes and store at -20°C for long-term storage. Avoid repeated freeze-thaw cycles.</li> </ol> |
| <b>Shipping</b>       | <p>Lyophilized Form: Shipped at ambient temperature (or with blue ice).</p> <p>Liquid Form: Shipped with blue ice (or on dry ice, depending on stability).</p> <p>Note: Upon receipt, please store the product immediately under the recommended conditions specified on the Certificate of Analysis (CoA).</p>                                                                                                                                                                                                                                                                                                                       |
| <b>Lead Time</b>      | <p>In-stock Items: Ships within 1-2 business days.</p> <p>Out-of-stock / Custom Orders: 2-4 weeks (Our team will confirm the exact timeline upon order placement).</p>                                                                                                                                                                                                                                                                                                                                                                                                                                                                |

This protein can be fully tailored to your specifications. Our platform directly supports expression host transitions, tag removal, and scalable bulk production.