



Recombinant pgi protein from *Lactobacillus casei*

Cat. No.: LBGF-0626-G992

This product is for research use only and is not intended for diagnostic use.

BASIC INFORMATION

Product Overview

Our Recombinant pgi protein is sourced from *Lactobacillus casei* and expressed via specialized heterologous hosts (*E. coli*, yeast, or mammalian cells) to achieve native-like conformation. We prioritize exceptional purity and low endotoxin levels to guarantee reliable performance in your experimental workflows. Detailed quality metrics, including purity and specific biological activity (where available), can be found on the lot-specific Certificate of Analysis (CoA).

Target	pgi protein
Source Organism	<i>Lactobacillus casei</i> ; strain BL23
Gene/Protein Name	pgi
Protein Sequence	MSYIKFDSSK LDKFVHANEL EQMQPLVTAA DKELREGTGA GKDFRGFIDL PVNYDKDEFA RIKAAAKKVQ GNSQVFVAIG IGGSYLGARM AVDFLSQTFR NLDPDLKFPE VYFAGNSISG TYLADLLDII GDRDFSINVI SKSGTTTEPS IAFRVLKAKL IEKYGKDGAK ERIYATTDRA KGALKQEADA EGYEEFVVPD DVGGRFSVMS AVGLLPPIAVA GGDIDEMMRG LGDGRKAYAS ADLKENEAYQ YAALRNILYR KGYTTTELLEN YEPTLQYLGE WWKQLMGESE GKDQKGIYPS

SANFSTDLHS LGQYIQEGLR NLMETVVWVE EPNRDLTIPE DANNLDGLGY LAGKKMSF
VN

RKAYEGVULA HTDGGVPVMT VSIPKQDAYT LGYLIYFFEA AVSISGYLNG INPFNQPGVE

AYKKNMFALL GRPGYEDMTK ELNARL

Subcellular Localization Cytoplasm

Protein Region Full Length

UniProt Accession Number B3WDC0

Biological Activity Not Test

Protein Family GPI family

Applications ELISA, Western Blot (WB), Immunoprecipitation (IP), Positive Control.
Note: Optimal dilutions/concentrations should be determined by the end user.

CHEMICAL & PHYSICAL PROPERTIES

Tag&Position Determined during the production process by default. If a specific tag is required, please notify us, and its development will be prioritized.

Purity ≥85% (SDS-PAGE)

Formulation

1. Liquid Form: Supplied in a default Tris/PBS-based storage buffer supplemented with 5% -50% glycerol.
2. Lyophilized Form: Lyophilized from a Tris/PBS-based pre-lyophilization buffer (pH 8.0) containing 6% Trehalose as a lyoprotectant.

Format Liquid or Lyophilized powder

HANDLING & STORAGE

Stability&Storage

1. Lyophilized Powder: Stored at -20°C to -80°C, the product is stable for 12 months from the date of receipt. Upon reconstitution, it is stable for up to 1 week at 2-8°C or 3 months at -20°C to -80°C.
2. Liquid Form: Stored at -20°C to -80°C, the product is stable for 6 to 12 months (avoid repeated freeze-thaw cycles). Aliquot upon arrival to prevent degradation.

Reconstitution	For Lyophilized Powder. Always centrifuge tubes before opening to bring the contents to the bottom. Do not vortex. 1. Reconstitute in sterile ddH₂O to a stock concentration of 0.1-1.0 mg/mL. 2. Dissolve the protein by gentle pipetting up and down or by inverting the tube. 3. Allow the solution to stand at room temperature for at least 20 minutes to ensure complete reconstitution. 4. Aliquot the reconstituted solution into working volumes and store at -20°C for long-term storage. Avoid repeated freeze-thaw cycles.
Shipping	Lyophilized Form: Shipped at ambient temperature (or with blue ice). Liquid Form: Shipped with blue ice (or on dry ice, depending on stability). Note: Upon receipt, please store the product immediately under the recommended conditions specified on the Certificate of Analysis (CoA).
Lead Time	In-stock Items: Ships within 1-2 business days. Out-of-stock / Custom Orders: 2-4 weeks (Our team will confirm the exact timeline upon order placement).

This protein can be fully tailored to your specifications. Our platform directly supports expression host transitions, tag removal, and scalable bulk production.