



Bacillus atrophaeus Spore Suspensions 9372, 10⁸

Cat. No.: LBGF-0926-GF19

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	<i>Bacillus atrophaeus</i> Spore Suspensions is a calibrated bacterial spore suspension suitable for direct inoculation of products or selection of substrates to monitor sterilization processes.
Target	<i>Bacillus atrophaeus</i>
Size	5 mL per vial
Population	10 ⁸ CFU/mL
Format	Suspensions
Diluent	Aqueous or 20% ethanol
Packaging	Penicillin bottle
Storage	2-8°C
Shelf Life	12 months from date of manufacture.
Applications	Used to validate sterilizers and monitor the effectiveness of the sterilization process.